Z-DNE® get back, faster





FAST RECOVERY COLD & AIR COMPRESSION THERAPY UNIT





NO ICE / NO WATER Therapeutic cooling technology without ice or water inside.

CLOSED & CLEAN SYSTEM

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off the liquid.

INNOVATIVE COOLING

Advance technology based on real-time temperature parameters.

STAY FOCUSED ON TEMPERATURE Set exact temperature you want to reach.

INTERMITTENT COMPRESSION

Improves benefits of cold therapy and healing with 3 levels of AIR compression.

NEW SYSTEM OF COMPRESSION MANAGEMENT Posibility to select different AIR compression setup.

INSPIRED BY YOUR SHAPE Complete range of Anatomic wraps designed on your body.

LIGHTER & STRONGER





small and powerful than ever

PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

FRIENDLY USE

All treatment parameters manageable in one screen.

ALWAYS ON THE MOVE

Light and handy can easily be with you anytime, anywhere.

REPEAT YOUR FAVORITE Possibility to store your personal treatments for frequent use.





Z-ONE MG455B







Designed and Delveloped by HTFenix®GmbH - Austria

ZAMAR MEDICAL d.o.o.

Sv. Martin, 6 - 52450 - Vrsar - (HR) Tel: +385 (0)52 496 111 - VAT: 96642432252 info@zamar.care

Features

APPLICATION	COLD
TEMPERATURE RANGE	+5°C +15°C (41°F 59°F)
SCREEN	LCD 5" color TouchScreen
WRAPS	Glycol Precharged
Glycol	NON TOX freeze
Air Compression levels	mmHg: 40 - 60 - 80
Working time	up to 10 h.
Medical class for Applied Parts	BF
Number of Wraps	up to 2
Storable therapies	1 + 2
DIMENSIONS	cm. 29 x 26 x 30h inch 11 x 10 x 12h
WEIGHT	9 kg (20 lb.)

- LCD 5", TFT color touch screen
- Easy Software, Intuitive Menu icons
- Possibility to select temperature and treatment time
- Possibility to manage different levels and Duration of the AIR compression
- 2 storable programs that can be customized by user
- Totally automated cycle management
- Fast connection of wraps and pipes with Dripless System
- Anatomic Wraps pre charge with glycol
- NON-TOX cooling glycol

0CATMG455BA00E EN Rev.00 2021.09

www.zamar.care





